



Assignment #5 - Recrystallization Worksheet

i You cannot edit this entry after it is graded.

Description Complete observations while in lab. Due at 5:00 pm the day following your lab section.

I worked in a group with

The work for this assignment is in My notebook

Grade 10 / 10

Graded on Sep 05, 2024, 12:32 PM CDT

Ilana Berlin - Sep 03, 2024, 2:31 PM CDT

Solvent Selection and Recrystallization of your impure solid

Name of your solid: 2-iodobenzoic acid

literature mp of your solid (°C): 160 - 164°C

Describe the original appearance of the solid:

opaque dull white powder

Table of Observations for Solubility Tests

solvent	bp of solvent, °C	rt	hot	re-cooled
water	100°C	ins	ins	ins
methanol	65°C	sol	N/A	N/A

ethanol	78°C	sol	N/A	N/A
ethyl acetate	77°C	sol	N/A	N/A
toluene	110.6C	ins	sol	ins
hexane	69°C	ins	ins	ins
solvent pair:	NA	ins	sol	ins

table abbreviations: sol = soluble; ins = insoluble; sl sol = slightly soluble (do not substitute)

Solvents that could be used: Methanol, Ethanol, Ethyl Acetate, Toluene

Solvent or solvent pair chosen and explain choice: Toluene, the only solvent that dissolved at BP instead of at room temperature. This will allow for recrystallization.

Mass of compound before recrystallization (mg): 110 mg

Approximate amount of solvents: 1.2 mL

Describe the appearance after recrystallization:

Small opaque dull white flakes that formed at the bottom of the tube and along the wooden stirring stick

Mass of compound after recrystallization (mg): 65.0

Percent Recovery: 59.0%

Observed mp (°C): 154.5 - 156.8 °C

Short Answer Questions:

Thinking back over both recrystallization part 1 and part 2 please complete the following short answer questions.

1. What is the purpose of recrystallization?

Recrystallization helps purify a compound.

2. At the end of recrystallization, where should the impurities be located?

The impurities should remain in the solvent.

3. Explain how the following items would adversely affect the percent recovery or purity of a recrystallization:

a. An unnecessarily large volume of solvent is used to dissolve the crude solid

Less of the solid will recrystallize out of the solution if too much solvent is used. This will decrease the percent recovery but will not effect the purity.

b. The crystals obtained after filtration are not washed with fresh cold solvent before drying

If the crystal are washed with solvent but the solvent is not cold them some of the crystals may redissolve, decreasing the percent recovery.

If the crystals are not washed at all then it may be harder to collect the crystals, decreasing the percent recovery.

c. The crystals obtained after filtration are washed with fresh hot solvent

The crystals will redissolve in the solvent and percent recovery will decrease.

d. Crystallization is accelerated by immediately placing the flask of hot solution in an ice bath

Cooling the sample too quickly then the impurities may also solidify, reducing the purity.

Erica Schultz - Aug 16, 2021, 1:12 PM CDT

Annotated IR spectra

Ilana Berlin - Sep 03, 2024, 2:50 PM CDT

