



## Assignment #12 - LabReport6

You cannot edit this entry after it is graded.

Description Beer's Law of Cu(II) ions in coordination complex.

I worked in a group with

The work for this assignment is in My notebook

**Grade** 9 / 10

Graded on Oct 06, 2023, 9:47 AM CDT

Jason Cody - Oct 01, 2021, 10:36 AM CDT

**TITLE:** (insert experimental title here. All italicized text in parentheses should be followed and then deleted throughout this template).

**Purpose:** (insert experimental purpose here).

**Reference:** Kotley, L. J., *Introduction to Chemistry in the Laboratory*, 20<sup>th</sup> Ed., Lake Forest College, 2023, Experiment xx, Appendix xx. (Edit the experiment title and/or appendix letter; add other references, if used, following the same format).

**Observation and Data:** (Write your color, concise, complete, past tense, passive voice description or narrative of the experiment as the experiment is performed. Complete sentences are used throughout.)

(If needed, insert tables and edit the header: Table 1. Preparation of Standard Solutions. If needed, insert figures and edit the caption below the figure: Figure 1. Beer's Law Plot of #12 Standard Solutions at  $\lambda = 520$  nm. Number tables and figures in order of appearance in the report.)

**Calculations:** (insert sample calculation here, if relevant. Otherwise, delete this section entirely).

**Conclusion:** (restate the quantitative values (percent error and/or CV) to indicate how well the goals of the experiment have been met; answer any questions in the experimental instructions, etc).

ReportTemplate.docx (15.5 kB)

Ilana Berlin - Oct 02, 2023, 11:29 AM CDT

## Date and Title

Ilana Berlin - Oct 02, 2023, 11:39 AM CDT

Purity and concentration of copper(II) ions in copper ammonium coordination compound.

## Purpose

Jason Cody - Oct 06, 2023, 9:42 AM CDT

By diluting known concentrations of copper ammonium coordination compound and using Beer's Law the concentration (M) of copper in a synthesized copper ammonium coordination compound **which complex, specifically? Did you make it last week?** can be determined. **what will you determine, exactly?**

Jason Cody - Sep 22, 2020, 3:19 PM CDT

## Reference

Ilana Berlin - Oct 02, 2023, 2:48 PM CDT

Kateley, L. J., *Introduction to Chemistry in the Laboratory, 20th Ed.*, Lake Forest College, 2021, Experiment 6, Appendix B,C.

Jason Cody - Sep 22, 2020, 3:19 PM CDT

## Data and Observations

Jason Cody - Oct 06, 2023, 9:45 AM CDT

Mass of dried blue powder copper(II) ammonium coordination compound, weighed using top loading balance, has a mass of 18.628g with the base of the petri dish (17.333g). The total weight of the synthesized compound is 1.295g.

10mL of a pale blue 5.4053g/L Cu(II) solution were added to a 30mL Pyrex beaker. 30mL of a 1M nitric acid solution was prepared using 5mL of a 6M nitric acid solution and 25mL of water. 15mL of water were added to a 50mL Pyrex beaker, then 5mL nitric acid were added, then the remaining 10mL of water. Five samples of varying concentration were prepared (Figure 1). Vial 1, containing just nitric acid solution, is completely clear. Vial 5, containing just  $\text{Cu}(\text{NO}_3)_2$  solution is a pale sky blue liquid. Each sample get progressively more saturated. All samples are transparent. **OK**

A piece of dark blue/purple hard solid synthesized  $[\text{Cu}(\text{NH}_3)_4]\text{SO}_4 \cdot \text{H}_2\text{O}$  was crushed using mortar and pestle. Crushing the compound brought out the purple color. 0.1464g (measure using an analytical balance) of the compound were added to a 10.00mL Kimax volumetric flask. 10mL of 1M nitric acid was added to the flask. The flask was lightly shaken to dissolve the compound. Approximately 3000 $\mu\text{m}$  of solution was added to vial 6. The solution in vial 6 is most similar in color to vials 3 and 4. **OK**


A microLAB spectrometer, model FS-528, was used to scan samples (Figure 2). The  $\lambda_{\text{max}}$  for Cu(II) is 770nm. The Beer's Law equation given by the samples is  $15.29682x - 0.09073$  **OK**

Figure 1. Preliminary Data **still preliminary?**

Vial Number	Volume of $\text{Cu}(\text{NO}_3)_2$ solution ( $\mu\text{L}$ )	Volume of $\text{HNO}_3$ ( $\mu\text{L}$ )	Concentration (M)	Absorbance (A)
1 (Blank)	0	3000	0	N/A

2	750	2250	0.0213	0.2108
3	1500	1500	0.04253	0.5928
4	2250	750	0.0628	0.8763
5	3000	0	0.08506	1.1952
6 (Unknown)	N/A	N/A	0.0621	0.8597

Figure 2. Spectrometer Data



Jason Cody - Sep 22, 2020, 3:19 PM CDT

## Calculations

Jason Cody - Oct 06, 2023, 9:46 AM CDT

$$1.295\text{g} (1000\text{mg}/1\text{g})(1\text{mmol}/245.7\text{mg}) = 5.271 \text{ mmol of } [\text{Cu}(\text{NH}_3)_4]\text{SO}_4 \cdot \text{H}_2\text{O}$$

$$\text{Percent yield} = 100 \times (5.271/6.00) = 87.9\% \text{ OK. Very reasonable.}$$

$$5.4053\text{g/L Cu(II)} (1 \text{ mol Cu(II)}/63.546\text{g Cu(II)}) = 0.085061 \text{ mol/L} = 8.5061 \times 10^{-2} \text{ M}$$

$$1\text{M}(30\text{mL}) \times 6\text{M}(\text{xmL}) = 5 \text{ mL}$$

$$30\text{mL} - 5\text{mL} = 25\text{mL}$$

$$8.5061 \times 10^{-2} \text{ M}(750\mu\text{m}) = \text{xM}(3000\mu\text{m}) = 2.1265 \times 10^{-2}\text{M}$$

$$8.5061 \times 10^{-2} \text{ M}(1500\mu\text{m}) = \text{xM}(3000\mu\text{m}) = 4.2531 \times 10^{-2}\text{M}$$

$$8.5061 \times 10^{-2} \text{ M}(2250\mu\text{m}) = \text{xM}(3000\mu\text{m}) = 6.2796 \times 10^{-2}\text{M}$$

$0.1464\text{g} (1\text{mol}/245.7\text{g})(1\text{ mol Cu}/1\text{ mol } [\text{Cu}(\text{NH}_3)_4]\text{SO}_4\cdot\text{H}_2\text{O}) = 5.958\times 10^{-4}\text{ mol of Cu}/0.01000\text{L} = 0.05958\text{M} = 5.958\times 10^{-2}\text{M}$

Percent error -  $100\times (0.0621-0.05958)/0.05958 = 4.4\%$  OK; comparison of calculated concentration from line equation?

Jason Cody - Sep 22, 2020, 3:19 PM CDT

## Conclusions

Jason Cody - Oct 06, 2023, 9:47 AM CDT

The compound was relatively pure since the concentration was close to the highest known concentration solution. ?? how does this indicate purity? There may be unused reactants that contribute to the diminished purity. There may have also been moisture remaining in the sample and  $\text{NH}_3$  lost in the  $\text{Cu}(\text{II})$  complex.

The b value (0.09073) in the Beer's law equation shows how far off the entire line is from the y intercept being 0. Since this number is close to zero, my plot was close to ideal.  $R^2$  is also a measure of the quality of your line (and, therefore, pipet technique). There may have been operative errors in measurement. OK, based on what?