

## Acid Equilibrium Constant for Methyl Red

### Purpose: Did you use your original solutions?

- Use Beer's Law to determine absorbance and solve for concentration of [HIn]
- Use Beer's Law to determine absorbance and solve for concentration of [In<sup>-</sup>]
- Determine pKa using Henderson – Hasselbach equation  $\text{pH} = \text{pKa} + \log\left(\frac{[\text{In}^-]}{[\text{HA}]}\right)$
- Check accuracy and precision of pKa using percent error and CV
- $\text{HIn} \rightleftharpoons \text{H}^+ + \text{In}^-$

### Reference:

- (1) Kateley, L. J., *Introduction to Chemistry in the Laboratory*, 20<sup>th</sup> Ed., Lake Forest College, 2021, Experiment 7\_Acid\_Equilibrium\_for\_Methyl\_Red Appendix C\_Spectroscopy.

### Prepare Solutions of HIn

- Pink transparent solution that darkened with concentration
- Relatively even amounts in each vial
- Blank of deionized water




TABLE 1: HIn/ Red SOLUTIONS, pH 2

Solution	Relative conc., C	Original HIn, $\mu\text{L}$	0.010 M HCl, $\mu\text{L}$	A at 520 nm
1, blank	0	0	<i>Provided</i>	0
2	0.400	1200	1800	0.388
3	0.600	1800	1200	0.604
4	0.800	2400	600	0.722
5	1.000	Enough to read in vial	0	0.786



### Calculations:

$$A_{525} = kC + \text{intercept}$$


$$A_{525} = 0.8131[\text{HIn}] + 0.0449$$

$$\text{Fit} = 0.9951$$

- The fit is good, intercept small

### Prep In- Solutions

- Saturated yellow transparent solution that darkened with concentration
- Even amounts of solution in each vial
- Same blank of deionized water from HnI solutions used



**TABLE 2: In<sup>-</sup>/ Yellow SOLUTIONS, pH 8**

Solution	Relative conc., C	Original In <sup>-</sup> , μL	0.010 M NaOAc, μL	A at 425 nm
6, blank	0	0	<i>Provided</i>	0
7	0.400	1200	1800	0.234
8	0.600	1800	1200	0.345
9	0.800	2400	600	0.528
10	1.000	Enough to read in vial	0	0.786

**Calculations:**

$$A_{425} = k_{425}[\text{In}^-] + \text{intercept}$$


$$A_{425} = 0.6188[\text{In}^-] - 0.0054$$

$$\text{Fit} = 0.9951$$

- The fit is good, intercept small, close to zero

**Prep Buffer Solution of  $\text{In}^-$** 

- 0.020M HOAc, 0.04M NaOAc, and standardized methyl red of varying concentrations
- Colors ranged from pale red (highest concentration of HOAc) to golden yellow (lowest concentration of HOAc)

**TABLE 3: BUFFERED SOLUTIONS**

Soln	Std methyl red, $\mu\text{L}$	0.02 M HOAc, $\mu\text{L}$	0.04 M NaOAc, $\mu\text{L}$	$\text{H}_2\text{O}$ , $\mu\text{L}$	pH	A at 520 nm	A at 425 nm	[HIn] (Red)	$[\text{In}^-]$ (Yellow)	Calculated $\text{pK}_a$
11, blank	0	500	500	1000	NA	0	0	NA	NA	NA
12	2000	5000	2500	500	4.83	0.199	0.158	0.190	0.264	4.78
13	2000	2500	2500	3000	5.12	0.158	0.153	0.139	0.256	4.90
14	2000	1000	2500	4500	5.55	0.139	0.139	0.116	0.282	5.17
15	2000	500	2500	5000	5.77	0.113	0.185	0.0838	0.308	5.32

**Sample Calculations for Vial 12:**

Calculation 1 : [HIn]

$$A_{525} = 0.8131[\text{HIn}] + 0.0449$$

$$0.199 = 0.813[\text{HIn}] + 0.0449$$

$$[\text{HIn}] = 0.190\text{M}$$

Calculation 2:  $[\text{In}^-]$

$$A_{425} = 0.6188[\text{In}^-] - 0.0054$$

$$0.158 = 0.619[\text{In}^-] - 0.0054$$

$$[\text{In}^-] = 0.263\text{M}$$

Calculation 3: pKa

$$\text{pH} = \text{pKa} + \log([\text{In}^-]/[\text{HIn}])$$

$$4.90 = \text{pKa} + \log(0.264/0.190)$$

$$\text{pKa} = 4.78$$

Calculation 4: Overall calculations

$$\text{pKa mean} = (4.78 + 4.90 + 5.17 + 5.32)/4 = 5.04$$

$$\text{Standard deviation} = 0.247 \text{ } 0.25$$

$$\text{CV} = 0.2547/5.04 \times 100 = 4.9\%$$

$$\text{Percent error} = 5.04 - 5.05/5.05 \times 100 = 0.2\%$$

- Good accuracy and bad precision